

Tungsten Syringe Shields for PET Radiopharmaceuticals



Mol-Image[®]
International Molecular Imaging

Headquarters:

Mongeri Building 19 Mayıs Mh.
Dr. Sevket Bey Sk No:5
34360 Sisli-Istanbul / Turkey
Phone : +90 212-231 0303
Fax : +90 212-219 8222

R&D:

Tubitak-MAM Teknoloji
Gelistirme Serbest Bolgesi
41470 Gebze-Kocaeli / Turkey
Phone : +90 262-644 4020
Fax : +90 262-644 9295

- 18 g / cc high density tungsten shield
- 9 mm thickness for maximum protection
- 2,1 mm lead equivalent viewing window
- Original cam lock mechanism for single handed locking and unlocking
- Available in two models for 2cc (2,5 cc) and 5 cc syringes

